

**In the Specification:**

Please amend the specification as follows. On page 11, please replace the paragraph starting on line 8 with the following paragraph:

For sequencing of the DNA inserts from phage clones reacting with anti-GD3 ganglioside mAbs, double stranded DNA is extracted by using commercially available kits (QIA prep. Spin, Miniprep kit, Qiagen Inc., Valencia, CA) from phage infected *E. coli* culture. The DNA sequencing can be performed with suitable primers (such as 5'-GCCAATAGTAGCACCAACGA-3' (SEQ ID NO:7), which is complementary to the synthetic pVIII gene and locates at more than 90 bp. downstream of the DNA inserts). DNA sequencing is carried out by standard techniques.